

User manual by product category

Caution:

These instructions for use refer to products that are only to be used by qualified professionals (i.e. dentists and/or oral and maxillofacial surgeons) unless expressly indicated otherwise (i.e. for dental technicians). All personnel must be specially trained and the products must be properly handled to avoid any risk of infection. When necessary, the qualified personnel must explain the hygiene rules that the patient must adopt and plan a regular outpatient control program.

Product identification

Equival® compatible prosthetic components manufactured by Phentagon Lab S.r.l. are medical devices that are intended to be used in the oral cavity. The functions of the prosthetic components are: reconditioning of the gums, anchoring to dental implants to support dental prostheses and the taking of dental impressions.

Healing screws

They are used when the implant to which they are fixed is removed and are left in place until the gums and the suturing around the screws have suitably healed. These components allow the soft tissues above the dental implant to heal so that the temporary and permanent prosthetic abutments can be inserted. The product is supplied decontaminated in order to reduce the microbiological load. They are available in various diameters and heights suitable for different gingival tissue thicknesses. Long-term, surgical-grade, single-use, non-sterile invasive devices supplied in individual or multi-packs.

Abutments and components for traditional screwed type prostheses or for the Multi-Unit technique

Screw-in and pass-through abutments complete with fastening screws are devices made of grade 5 titanium which anchor the prosthesis by means of screws that pass through the dental implant. The product is supplied decontaminated in order to reduce the microbiological load. They are available in various diameters and heights suitable for different gingival tissue thicknesses. Long-term, surgical-grade, single-use, non-sterile invasive devices complete with fastening screws supplied in individual or multi-packs.

Temporary abutments and components for traditional screwed type prostheses or for the Multi-Unit technique

Temporary abutments complete with fastening screws are devices made of grade 5 titanium which anchor the prosthesis by means of screws that pass through the dental implant or using screws that are fixed to the Multi-Unit abutment. The product is supplied decontaminated in order to reduce the microbiological load. They are available in various diameters and heights to adapt to various clinical needs. Long-term, surgical-grade, single-use, non-sterile invasive devices supplied in individual or multi-packs.

Screws for Abutments, Temporary Abutments or for the Multi-Unit technique

These components are made of grade 5 titanium and are used to fix the prosthesis to the permanent implant. They are available in various diameters and heights to ensure compatibility with the dental implant. They are non-surgical grade, invasive, medical device accessories intended for long-term use in the oral cavity and are supplied with the respective abutments or sold in packages as spare parts. Single-use, non-sterile supplied in individual or multi-packs.

Sterilisation

The devices do not have an expiry date and are supplied NON-STERILE. The user must autoclave sterilise the devices using a validated cycle in compliance with the requirements of standard UNI EN ISO 17665-1:2007 (134°C, 10 min., 2 bar) before each use.

Storage

After sterilisation, the Equival® prosthetic components must be stored in the sterilisation pouches and only be opened immediately before use. Unless they become damaged, a sterile environment is normally maintained within the autoclave sterilisation pouches. Therefore, be sure not to use the components if the pouches containing the devices are damaged and to re-sterilise them in new pouches before use. The storage period for sterilised products inside the pouches must not exceed that recommended by the pouch manufacturer. The product must be stored in a cool and dry place, away from direct sunlight, water and sources of heat.

Warranty and limitations

The success of all implants requires that the surgical technique be meticulously performed. Careful patient selection and a final prosthetic restoration that meets the patient's individual needs are essential. The surgical planning of single or multiple implants and the selecting of the suitable shape and size of the prosthetic components for the patient's anatomy are also crucial for the success of the implant. Phentagon Lab s.r.l. therefore strongly recommends that all dentists take one or more endosseous dental implant courses to properly learn about the implant system. Since the factors relating to the services provided by the operator, the selection of the patient and the surgical and rehabilitative technique are out of its control, Phentagon Lab s.r.l. will not replace the product due to failures or other undesired reactions resulting from operator errors or the lack of operator information on the use of an implant, and cannot be held liable for any unsuitable conduct by the patient.








The prosthetic components are made of a medical grade titanium alloy. Numerous animal and human studies conducted over several years by researchers and reported in many scientific publications show that titanium is highly biocompatible. However, the prevention of fractures of an implant or injuries derived from trauma or other overloading such as, for example, a severe traumatic occlusion, cannot be guaranteed.

For these reasons Phentagon Lab s.r.l. can only provide guarantees relating to the correct use of its products.

This warranty replaces all other warranties, expressed or implied, including the marketability and adaptability to particular situations. Phentagon Lab s.r.l. only and exclusively replaces products that have been returned to Phentagon Lab s.r.l., have undergone the necessary inspections and are proven to be defective. Phentagon Lab s.r.l. is not liable for damages deriving from accidents or events caused by third parties and if the relative claim was unfounded or incorrect.

Unless expressly specified otherwise, the device is designed for single use only, and any attempt to recondition it for subsequent reuse compromises the quality of the materials and the effectiveness of the device. Phentagon Lab s.r.l. expressly prohibits the re-use of its single-use devices.

Package contents: 1 Product shown on the label, 1 User manual

Symbol	Description
	Caution, read the instructions for use
	Lot number
REF	Product identification code
	Non sterile. The user must autoclave sterilise the product using a validated cycle in compliance with standard UNI EN ISO 17665-1:2007
	Manufacturer
	Read the instructions for use
	Single-use device (where applied; see label)
	CE conformity marking



PHENTAGON LAB s.r.l. Via F.M. Malfatti, 75 02100 - Rieti – Italy Tel 0746 247025 Fax 0746 247025 – info@phentagonlab.it